

## HTT FIBER Y LIGHT GUIDE CABLE

for Xenon light source; w/o connectors

Y-Typ 2300mm; 3.5/4.8; 1x Proj / 2x Endo

tube colour lilac

**HTT**  
HIGH TEMPERATURE TECHNOLOGY LINE



**SKU:** 540-816-23LL

**Categories:** [Endoscopy](#), [Equipment](#), [Laparoscopy](#), [Gynecology](#), [GUBBINI/Gyn](#), [Urology](#), [Arthroscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Light Guide Cables](#)

**TYPE** [HTT \(High Temperature Technology\)](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Light Guide Cables](#)

**TYPE** [HTT \(High Temperature Technology\)](#)