

## HTT FIBER Y LIGHT GUIDE CABLE

for Xenon light source; w/o connectors

Y-Typ 1800mm; 3.5/4.8; 2x Proj / 1x Endo

tube colour grey

HTT

HIGH TEMPERATURE TECHNOLOGY LINE



**SKU:** 540-817-18

**Categories:** [Endoscopy](#), [Equipment](#), [Laparoscopy](#), [Gynecology](#), [GUBBINI/Gyn](#), [Urology](#), [Arthroscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Light Guide Cables</a>
<b>TYPE</b>	<a href="#">HTT (High Temperature Technology)</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Light Guide Cables</a>
<b>TYPE</b>	<a href="#">HTT (High Temperature Technology)</a>