

## MINI-ENDOSCOPE HD+

Ø 1.9mm; 0°; 60 mm

autoclavable



**SKU:** 601-122-00HD

**Categories:** [Endoscopy](#), [Equipment](#), [Arthroscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

|              |   |
|--------------|---|
| <b>UNIT</b>  | <a href="#">pcs.</a>                              |
| <b>RANGE</b> | <a href="#">Telescopes</a>                        |
| <b>SCOPE</b> | <a href="#">Arthroscopy</a> , <a href="#">ENT</a> |



### ZUSÄTZLICHE INFORMATIONEN

|              |   |
|--------------|---|
| <b>UNIT</b>  | <a href="#">pcs.</a>                              |
| <b>RANGE</b> | <a href="#">Telescopes</a>                        |
| <b>SCOPE</b> | <a href="#">Arthroscopy</a> , <a href="#">ENT</a> |