

## HYSTEROSCOPE HD+

Ø 4.0mm; 0°; 300 mm

autoclavable

**SKU:** 601-114-00HD

**Categories:** [Endoscopy](#), [Equipment](#), [Gynecology](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Telescopes](#)

**SCOPE** [Gynecology](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Telescopes](#)

**SCOPE** [Gynecology](#)