

HYSTEROSCOPE HD+2

Ø 4.0mm; 30°; 300 mm

autoclavable

SKU: 601-114-30HD2

Categories: [Endoscopy](#), [Equipment](#), [Gynecology](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Telescopes](#)

SCOPE [Gynecology](#)



30°

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Telescopes](#)

SCOPE [Gynecology](#)