

HYSTEROSCOPE HD+2

Ø 4.0mm; 70°; 300 mm

autoclavable

SKU: 601-114-70HD2

Categories: [Endoscopy](#), [Equipment](#), [Gynecology](#)




ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Telescopes](#)

SCOPE [Gynecology](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Telescopes](#)

SCOPE [Gynecology](#)