

## HYSTEROSCOPY DIAGNOSTIC SHEATH

30°; 1 turnable stopcock

Ø 2.9 mm; 228 mm



**SKU:** 640-010-23

**Categories:** [Endoscopy](#), [Gynecology](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Hysteroscopy</a>
<b>TYPE</b>	<a href="#">Diagnostic Sheaths</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Hysteroscopy</a>
<b>TYPE</b>	<a href="#">Diagnostic Sheaths</a>