

## HYSTEROSCOPY DIAGNOSTIC SHEATH

30°; 1 turnable stopcock

Ø 4.0 mm; 270 mm; white

for telescope 2,9mm/30°/303mm



**SKU:** 640-010-33

**Categories:** [Endoscopy](#), [Gynecology](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Hysteroscopy</a>
<b>TYPE</b>	<a href="#">Diagnostic Sheaths</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Hysteroscopy</a>
<b>TYPE</b>	<a href="#">Diagnostic Sheaths</a>