

## INSTRUMENT BOX WITH LID WITH KNOB

300 x 200 x 50mm



**SKU:** 250-260-30

**Categories:** [Open Surgery](#), [Various Items](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bowls</a>
<b>KLS Martin</b>	56-411-30

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bowls</a>
<b>KLS Martin</b>	56-411-30