

IRRIGATION ELECTRODE; MONOPOLAR

trumpet valve; Luer-Lock

Ø 5.0mm; 330mm



SKU: 622-030-53

Categories: [Endoscopy](#), [Laparoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [HF-Electrodes, monopolar](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [HF-Electrodes, monopolar](#)