

## IRRIGATION ELECTRODE; MONOPOLAR

trumpet valve; Luer-Lock

Ø 5.0mm; 330mm



**SKU:** 622-030-53

**Categories:** [Endoscopy](#), [Laparoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [HF-Electrodes, monopolar](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [HF-Electrodes, monopolar](#)