

J-HOOK ELECTRODE; ADAPTABLE

for pistol handle 620-010-00

Ø 5.0mm; 330mm



SKU: 622-050-23

Categories: [Endoscopy](#), [Laparoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [HF-Electrodes, monopolar](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [HF-Electrodes, monopolar](#)