

JAMISON NEEDLE HOLDER

straight; serrated

23.5cm (9.4")



SKU: 050-148-23

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM270
KLS Martin	20-404-23
AUTHORS	JAMISON
Medicon	10.20.23
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM270
KLS Martin	20-404-23

AUTHORS

[JAMISON](#)

Medicon

10.20.23

SPECIFICATIONS

[serrated](#), [straight](#)