


## JAMISON NEEDLE HOLDER

straight; serrated

23.5cm (9.4")



**SKU:** 050-148-23

**Categories:** [Open Surgery](#), [Suture](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>Aesculap</b>	BM270
<b>KLS Martin</b>	20-404-23
<b>AUTHORS</b>	<a href="#">JAMISON</a>
<b>Medicon</b>	10.20.23
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>Aesculap</b>	BM270
<b>KLS Martin</b>	20-404-23

**AUTHORS**

[JAMISON](#)

**Medicon**

10.20.23

**SPECIFICATIONS**

[serrated](#), [straight](#)