

## JANET WOUND AND BLADDER SYRINGE

with two exchangeable tips; 200 ml

**SKU:** 087-070-20

**Categories:** [Open Surgery](#), [Puncture / Special Needles](#)



### ZUSÄTZLICHE INFORMATIONEN

|                   |                          |
|-------------------|--------------------------|
| <b>UNIT</b>       | <a href="#">pcs.</a>     |
| <b>TYPE</b>       | <a href="#">Syringes</a> |
| <b>Aesculap</b>   | SH420                    |
| <b>KLS Martin</b> | 18-930-20                |
| <b>AUTHORS</b>    | <a href="#">JANET</a>    |
| <b>Medicon</b>    | 09.78.85                 |

### ZUSÄTZLICHE INFORMATIONEN

|                   |                          |
|-------------------|--------------------------|
| <b>UNIT</b>       | <a href="#">pcs.</a>     |
| <b>TYPE</b>       | <a href="#">Syringes</a> |
| <b>Aesculap</b>   | SH420                    |
| <b>KLS Martin</b> | 18-930-20                |
| <b>AUTHORS</b>    | <a href="#">JANET</a>    |
| <b>Medicon</b>    | 09.78.85                 |