

JANSEN EAR FORCEPS

bayonet shaped; serrated

22.0cm (8.8"); width of jaw 2.3 +/-0.2



SKU: 020-760-22
 Categories: [Open Surgery](#), [Forceps](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Ear and Nose Dressing Forceps
AUTHORS	GRUENWALD
SPECIFICATIONS	bajonett , serrated



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Ear and Nose Dressing Forceps
AUTHORS	GRUENWALD
SPECIFICATIONS	bajonett , serrated