

JANSEN SELFRETAINING RETRACTOR

blunt; 3X3 teeth

10.0cm (4")


SKU: 070-631-10

Categories: [Open Surgery](#), [Retractors](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
Aesculap	BV032
KLS Martin	15-723-10
AUTHORS	JANSEN
Medicon	24.13.20
SPECIFICATIONS	blunt , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
Aesculap	BV032
KLS Martin	15-723-10

AUTHORS

[JANSEN](#)

Medicon

24.13.20

SPECIFICATIONS

[blunt](#), [toothed](#)