

JEWELER'S BIPOLAR COAGULATION FORCEPS

straight; pointed; Ø 0.25mm

11.0cm (4.4"); flat plug



SKU: 022-220-25
Categories: [Open Surgery](#), [Forceps](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
Medicon	88.12.05
SPECIFICATIONS	sharp , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
Medicon	88.12.05
SPECIFICATIONS	sharp , straight