

JOHNS-HOPKINS BULLDOG CLAMP

straight; length: 63mm

SKU: 026-018-06

Categories: [Hemostats](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Bulldog Clamps](#)

KLS Martin 13-119-06

AUTHORS [HOPKINS, JOHNS](#)

Medicon 15.15.63

SPECIFICATIONS [straight](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Bulldog Clamps](#)

KLS Martin 13-119-06

AUTHORS [HOPKINS, JOHNS](#)

Medicon 15.15.63

SPECIFICATIONS [straight](#)