

JOHNS-HOPKINS BULLDOG CLAMP

straight; length: 90mm

SKU: 026-018-09

Categories: [Open Surgery](#), [Hemostats](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Bulldog Clamps](#)

Aesculap BH017

AUTHORS [HOPKINS, JOHNS](#)

SPECIFICATIONS [straight](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Bulldog Clamps](#)

Aesculap BH017

AUTHORS [HOPKINS, JOHNS](#)

SPECIFICATIONS [straight](#)