

JOHNS-HOPKINS BULLDOG CLAMP

curved; length: 38mm

SKU: 026-019-04

Categories: [Open Surgery](#), [Hemostats](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
Aesculap	BH011
AUTHORS	HOPKINS, JOHNS
Medicon	15.15.38
SPECIFICATIONS	curved



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
Aesculap	BH011
AUTHORS	HOPKINS, JOHNS
Medicon	15.15.38
SPECIFICATIONS	curved