

JOHNS HOPKINS BULLDOG CLAMP

straight; length: 38mm

SKU: 026-018-04

Categories: [Open Surgery](#), [Hemostats](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
Aesculap	BH011
AUTHORS	HOPKINS, JOHNS
SPECIFICATIONS	straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
Aesculap	BH011
AUTHORS	HOPKINS, JOHNS
SPECIFICATIONS	straight