

## JOHNS-HOPKINS HAEMOSTATIC CLAMP

curved; serrated

20.0cm (8")




**SKU:** 040-497-20

**Categories:** [Open Surgery](#), [Hemostats](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Artery Forceps</a> , <a href="#">Dissecting and Ligature Forceps</a>
<b>KLS Martin</b>	13-467-20
<b>AUTHORS</b>	<a href="#">HOPKINS</a> , <a href="#">JOHNS</a>
<b>Medicon</b>	17.07.45
<b>SPECIFICATIONS</b>	<a href="#">curved</a> , <a href="#">serrated</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Artery Forceps</a> , <a href="#">Dissecting and Ligature Forceps</a>
<b>KLS Martin</b>	13-467-20
<b>AUTHORS</b>	<a href="#">HOPKINS</a> , <a href="#">JOHNS</a>
<b>Medicon</b>	17.07.45

## SPECIFICATIONS

[curved](#), [serrated](#)