

JOHNSON NEEDLE HOLDER

curved; serrated; grooved jaw

27.0cm (10.8")


SKU: 050-193-27

Categories: [Open Surgery](#), [Suture](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
AUTHORS	JOHNSON
Medicon	10.43.87
SPECIFICATIONS	curved , serrated



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
AUTHORS	JOHNSON
Medicon	10.43.87
SPECIFICATIONS	curved , serrated