

JONES TOWEL FORCEPS

sharp

5.0cm (2")



SKU: 030-105-05

Categories: [Open Surgery](#), [Hemostats](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Towel Forceps
Aesculap	BF401
KLS Martin	13-881-05
AUTHORS	JONES
Medicon	17.50.05
SPECIFICATIONS	curved , sharp

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Towel Forceps
Aesculap	BF401
KLS Martin	13-881-05

AUTHORS

[JONES](#)

Medicon

17.50.05

SPECIFICATIONS

[curved](#), [sharp](#)