

KERRISON MODIFIZ. PUNCH; 40° UP CUT

3mm; without ejector; thin footplate

20.0cm (8"); 15mm open




SKU: 240-272-03T

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
AUTHORS	KERRISON
SPECIFICATIONS	40° up cut , without ejector



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
AUTHORS	KERRISON
SPECIFICATIONS	40° up cut , without ejector