

SWING-SYSTEM ® KERRISON PUNCH

40° up cut; 5mm; without ejector

18.0cm (7.2"); 10mm open; CC




SKU: 249-071-05CC

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
MATERIAL	Ceramic(CC)
AUTHORS	KERRISON
Medicon	57.17.64
SPECIFICATIONS	40° up cut , without ejector



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
MATERIAL	Ceramic(CC)
AUTHORS	KERRISON
Medicon	57.17.64

SPECIFICATIONS

[40° up cut, without ejector](#)