

SWING-SYSTEM ® KERRISON PUNCH

40° up cut; 5mm; without ejector

thin foot; 18.0cm (7.2"); 10mm open; CC



SKU: 249-071-05TCC

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
MATERIAL	Ceramic(CC)
AUTHORS	KERRISON
Medicon	57.17.79
SPECIFICATIONS	40° up cut , thin footplate , without ejector



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
MATERIAL	Ceramic(CC)
AUTHORS	KERRISON
Medicon	57.17.79

SPECIFICATIONS

[40° up cut](#), [thin footplate](#), [without ejector](#)