

## SWING-SYSTEM ® KERRISON PUNCH

40° up cut; 1mm; without ejector

20.0cm (8"); 10mm open; CC



**SKU:** 249-072-01CC

**Categories:** [Open Surgery](#), [Neurosurgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Laminectomy Punches</a>
<b>MATERIAL</b>	<a href="#">Ceramic(CC)</a>
<b>AUTHORS</b>	<a href="#">KERRISON</a>
<b>Medicon</b>	57.17.68
<b>SPECIFICATIONS</b>	<a href="#">40° up cut</a> , <a href="#">without ejector</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Laminectomy Punches</a>
<b>MATERIAL</b>	<a href="#">Ceramic(CC)</a>
<b>AUTHORS</b>	<a href="#">KERRISON</a>
<b>Medicon</b>	57.17.68

## SPECIFICATIONS

[40° up cut, without ejector](#)