

SWING-SYSTEM® KERRISON PUNCH

40° up cut; 3mm; with ejector

20.0cm (8"); 10mm open; CC




SKU: 249-072-03ACC

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
Aesculap	FK914R
MATERIAL	Ceramic(CC)
AUTHORS	KERRISON
SPECIFICATIONS	40° up cut , with ejector



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
Aesculap	FK914R
MATERIAL	Ceramic(CC)
AUTHORS	KERRISON

SPECIFICATIONS

[40° up cut, with ejector](#)