

SWING-SYSTEM ® KERRISON PUNCH

40° up cut; 3mm; without ejector

thin foot; 20.0cm (8"); 10mm open; CC



SKU: 249-072-03TCC

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
MATERIAL	Ceramic(CC)
AUTHORS	KERRISON
Medicon	57.17.82
SPECIFICATIONS	40° up cut , thin footplate , without ejector



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
MATERIAL	Ceramic(CC)
AUTHORS	KERRISON
Medicon	57.17.82

SPECIFICATIONS

[40° up cut](#), [thin footplate](#), [without ejector](#)