

SWING-SYSTEM ® KERRISON PUNCH

40° up cut; 5mm; without ejector

20.0cm (8"); 10mm open; CC



SKU: 249-072-05CC

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
MATERIAL	Ceramic(CC)
AUTHORS	KERRISON
Medicon	249-072-05CC
SPECIFICATIONS	40° up cut , without ejector



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
MATERIAL	Ceramic(CC)
AUTHORS	KERRISON
Medicon	249-072-05CC

SPECIFICATIONS

[40° up cut, without ejector](#)