

SWING-SYSTEM ® KERRISON PUNCH

40° up cut; 2mm; without ejector

thin foot; 23.0cm (9.2"); 10mm open



SKU: 249-073-02T

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN


UNIT [pcs.](#)

TYPE [Laminectomy Punches](#)

AUTHORS [KERRISON](#)

Medicon 57.17.37

SPECIFICATIONS [40° up cut](#), [thin footplate](#), [without ejector](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Laminectomy Punches](#)

AUTHORS [KERRISON](#)

Medicon 57.17.37

SPECIFICATIONS [40° up cut](#), [thin footplate](#), [without ejector](#)