

SWING-SYSTEM ® KERRISON PUNCH


40° up cut; 2mm; without ejector

thin foot; 23.0cm (9.2"); 10mm open; CC



SKU: 249-073-02TCC

Categories: [Neurosurgery](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
MATERIAL	Ceramic(CC)
AUTHORS	KERRISON
Medicon	999.40.68
SPECIFICATIONS	40° up cut , thin footplate , without ejector

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
MATERIAL	Ceramic(CC)
AUTHORS	KERRISON
Medicon	999.40.68

SPECIFICATIONS

[40° up cut](#), [thin footplate](#), [without ejector](#)