

SWING-SYSTEM ® KERRISON PUNCH


40° up cut; 2mm; without ejector

thin foot; 23.0cm (9.2"); 10mm open; CC



SKU: 249-073-02TCC

Categories: [Open Surgery](#), [Neurosurgery](#)



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Laminectomy Punches |
| MATERIAL | Ceramic(CC) |
| AUTHORS | KERRISON |
| Medicon | 999.40.68 |
| SPECIFICATIONS | 40° up cut , thin footplate , without ejector |

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------|-------------------------------------|
| UNIT | pcs. |
| TYPE | Laminectomy Punches |
| MATERIAL | Ceramic(CC) |
| AUTHORS | KERRISON |
| Medicon | 999.40.68 |

SPECIFICATIONS

[40° up cut](#), [thin footplate](#), [without ejector](#)