

## SWING-SYSTEM ® KERRISON PUNCH

40° up cut; 3mm; with ejector; thin foot

23.0cm (9.2"); 10mm open; CC



**SKU:** 249-073-03ATCC

**Categories:** [Open Surgery](#), [Neurosurgery](#)

### ZUSÄTZLICHE INFORMATIONEN

|                       |  |
|-----------------------|--|
| <b>UNIT</b>           | <a href="#">pcs.</a>   |
| <b>TYPE</b>           | <a href="#">Laminectomy Punches</a>  |
| <b>MATERIAL</b>       | <a href="#">Ceramic(CC)</a>  |
| <b>AUTHORS</b>        | <a href="#">KERRISON</a>   |
| <b>SPECIFICATIONS</b> | <a href="#">40° up cut</a> , <a href="#">thin footplate</a> , <a href="#">with ejector</a> |



### ZUSÄTZLICHE INFORMATIONEN

|                       |  |
|-----------------------|--|
| <b>UNIT</b>           | <a href="#">pcs.</a>   |
| <b>TYPE</b>           | <a href="#">Laminectomy Punches</a>  |
| <b>MATERIAL</b>       | <a href="#">Ceramic(CC)</a>  |
| <b>AUTHORS</b>        | <a href="#">KERRISON</a>   |
| <b>SPECIFICATIONS</b> | <a href="#">40° up cut</a> , <a href="#">thin footplate</a> , <a href="#">with ejector</a> |