

## SWING-SYSTEM ® KERRISON PUNCH


40° up cut; 3mm; without ejector

thin foot; 23.0cm (9.2"); 10mm open; CC



**SKU:** 249-073-03TCC

**Categories:** [Open Surgery](#), [Neurosurgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Laminectomy Punches</a>
<b>MATERIAL</b>	<a href="#">Ceramic(CC)</a>
<b>AUTHORS</b>	<a href="#">KERRISON</a>
<b>Medicon</b>	999.40.69
<b>SPECIFICATIONS</b>	<a href="#">40° up cut</a> , <a href="#">thin footplate</a> , <a href="#">without ejector</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Laminectomy Punches</a>
<b>MATERIAL</b>	<a href="#">Ceramic(CC)</a>
<b>AUTHORS</b>	<a href="#">KERRISON</a>
<b>Medicon</b>	999.40.69

## SPECIFICATIONS

[40° up cut](#), [thin footplate](#), [without ejector](#)