

## SWING-SYSTEM ® KERRISON PUNCH

40° up uct; 4mm; without ejector

23.0cm (9.2"); 10mm open



**SKU:** 249-073-04

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

|                       |  |
|-----------------------|--|
| <b>UNIT</b>           | <a href="#">pcs.</a>   |
| <b>TYPE</b>           | <a href="#">Laminectomy Punches</a>                          |
| <b>AUTHORS</b>        | <a href="#">KERRISON</a>                                     |
| <b>Medicon</b>        | 57.17.34   |
| <b>SPECIFICATIONS</b> | <a href="#">40° up cut</a> , <a href="#">without ejector</a> |



### ZUSÄTZLICHE INFORMATIONEN

|                       |  |
|-----------------------|--|
| <b>UNIT</b>           | <a href="#">pcs.</a>   |
| <b>TYPE</b>           | <a href="#">Laminectomy Punches</a>                          |
| <b>AUTHORS</b>        | <a href="#">KERRISON</a>                                     |
| <b>Medicon</b>        | 57.17.34   |
| <b>SPECIFICATIONS</b> | <a href="#">40° up cut</a> , <a href="#">without ejector</a> |