

## SWING-SYSTEM ® KERRISON PUNCH

90° up cut; 3mm; with ejector

20.0cm (8"); 10mm open; thin foot, CC



**SKU:** 249-052-03TCC

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Laminectomy Punches</a>
<b>MATERIAL</b>	<a href="#">Ceramic(CC)</a>
<b>AUTHORS</b>	<a href="#">KERRISON</a>
<b>SPECIFICATIONS</b>	<a href="#">90° up cut</a> , <a href="#">thin footplate</a> , <a href="#">with ejector</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Laminectomy Punches</a>
<b>MATERIAL</b>	<a href="#">Ceramic(CC)</a>
<b>AUTHORS</b>	<a href="#">KERRISON</a>
<b>SPECIFICATIONS</b>	<a href="#">90° up cut</a> , <a href="#">thin footplate</a> , <a href="#">with ejector</a>