

## SWING-SYSTEM® KERRISON PUNCH

40° up cut; 2mm; without ejector

18.0cm (7.2"); 10mm open; thin foot



**SKU:** 249-071-02T

**Categories:** [Open Surgery](#), [Neurosurgery](#)

### ZUSÄTZLICHE INFORMATIONEN

|                       |   |
|-----------------------|---|
| <b>UNIT</b>           | <a href="#">pcs.</a>  |
| <b>TYPE</b>           | <a href="#">Laminectomy Punches</a>   |
| <b>Aesculap</b>       | FG742   |
| <b>AUTHORS</b>        | <a href="#">KERRISON</a>  |
| <b>Medicon</b>        | 57.17.17  |
| <b>SPECIFICATIONS</b> | <a href="#">40° up cut</a> , <a href="#">thin footplate</a> , <a href="#">without ejector</a> |



### ZUSÄTZLICHE INFORMATIONEN

|                 |                                     |
|-----------------|-------------------------------------|
| <b>UNIT</b>     | <a href="#">pcs.</a>                |
| <b>TYPE</b>     | <a href="#">Laminectomy Punches</a> |
| <b>Aesculap</b> | FG742                               |
| <b>AUTHORS</b>  | <a href="#">KERRISON</a>            |
| <b>Medicon</b>  | 57.17.17                            |

## SPECIFICATIONS

[40° up cut](#), [thin footplate](#), [without ejector](#)