

SWING-SYSTEM® KERRISON PUNCH

40° up cut; 3mm; with ejector

18.0cm (7.2"); 10mm open




SKU: 249-071-03A

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Laminectomy Punches |
| Aesculap | FG823 |
| AUTHORS | KERRISON |
| Medicon | 57.17.13 |
| SPECIFICATIONS | 40° up cut , with ejector |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------|-------------------------------------|
| UNIT | pcs. |
| TYPE | Laminectomy Punches |
| Aesculap | FG823 |
| AUTHORS | KERRISON |
| Medicon | 57.17.13 |

SPECIFICATIONS

[40° up cut, with ejector](#)