

## SWING-SYSTEM ® KERRISON PUNCH

40° up cut; 3mm; with ejector; thin foot

18.0cm (7.2"); 10mm open; CC



**SKU:** 249-071-03ATCC

**Categories:** [Open Surgery](#), [Neurosurgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Laminectomy Punches</a>
<b>MATERIAL</b>	<a href="#">Ceramic(CC)</a>
<b>AUTHORS</b>	<a href="#">KERRISON</a>
<b>SPECIFICATIONS</b>	<a href="#">40° up cut</a> , <a href="#">thin footplate</a> , <a href="#">with ejector</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Laminectomy Punches</a>
<b>MATERIAL</b>	<a href="#">Ceramic(CC)</a>
<b>AUTHORS</b>	<a href="#">KERRISON</a>
<b>SPECIFICATIONS</b>	<a href="#">40° up cut</a> , <a href="#">thin footplate</a> , <a href="#">with ejector</a>