

SWING-SYSTEM ® KERRISON PUNCH

40° up cut; 3mm; without ejector

thin foot; 18.0cm (7.2"); 10mm open



SKU: 249-071-03T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN


UNIT [pcs.](#)

TYPE [Laminectomy Punches](#)

AUTHORS [KERRISON](#)

Medicon 57.17.18

SPECIFICATIONS [40° up cut](#), [thin footplate](#), [without ejector](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Laminectomy Punches](#)

AUTHORS [KERRISON](#)

Medicon 57.17.18

SPECIFICATIONS [40° up cut](#), [thin footplate](#), [without ejector](#)