

SWING-SYSTEM ® KERRISON PUNCH


40° up cut; 4mm; without ejector

18.0cm (7.2"); 10mm open; CC



SKU: 249-071-04CC

Categories: [Open Surgery](#), [Neurosurgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
MATERIAL	Ceramic(CC)
AUTHORS	KERRISON
Medicon	57.17.63
SPECIFICATIONS	40° up cut , without ejector

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
MATERIAL	Ceramic(CC)
AUTHORS	KERRISON
Medicon	57.17.63

SPECIFICATIONS

[40° up cut, without ejector](#)