

## KERRISON PUNCH; 40° UP CUT

3mm; with ejector; thin footplate

23.0cm (9.2"); 10mm open




**SKU:** 240-073-03AT

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

|                       |   |
|-----------------------|---|
| <b>UNIT</b>           | <a href="#">pcs.</a>                                      |
| <b>TYPE</b>           | <a href="#">Laminectomy Punches</a>                       |
| <b>AUTHORS</b>        | <a href="#">KERRISON</a>                                  |
| <b>SPECIFICATIONS</b> | <a href="#">40° up cut</a> , <a href="#">with ejector</a> |



### ZUSÄTZLICHE INFORMATIONEN

|                       |   |
|-----------------------|---|
| <b>UNIT</b>           | <a href="#">pcs.</a>                                      |
| <b>TYPE</b>           | <a href="#">Laminectomy Punches</a>                       |
| <b>AUTHORS</b>        | <a href="#">KERRISON</a>                                  |
| <b>SPECIFICATIONS</b> | <a href="#">40° up cut</a> , <a href="#">with ejector</a> |