

**KERRISON PUNCH; 40° UP CUT**

5mm; without ejector

23.0cm (9.2"); 10mm open




**SKU:** 240-073-05

**Categories:** [Neurosurgery](#), [Open Surgery](#)

**ZUSÄTZLICHE INFORMATIONEN**

|                       |                                                              |
|-----------------------|--------------------------------------------------------------|
| <b>UNIT</b>           | <a href="#">pcs.</a>                                         |
| <b>TYPE</b>           | <a href="#">Laminectomy Punches</a>                          |
| <b>Aesculap</b>       | FF685                                                        |
| <b>AUTHORS</b>        | <a href="#">KERRISON</a>                                     |
| <b>SPECIFICATIONS</b> | <a href="#">40° up cut</a> , <a href="#">without ejector</a> |



**ZUSÄTZLICHE INFORMATIONEN**

|                       |                                                              |
|-----------------------|--------------------------------------------------------------|
| <b>UNIT</b>           | <a href="#">pcs.</a>                                         |
| <b>TYPE</b>           | <a href="#">Laminectomy Punches</a>                          |
| <b>Aesculap</b>       | FF685                                                        |
| <b>AUTHORS</b>        | <a href="#">KERRISON</a>                                     |
| <b>SPECIFICATIONS</b> | <a href="#">40° up cut</a> , <a href="#">without ejector</a> |