

## KERRISON PUNCH; 40° UP CUT

3mm; without ejector

18.0cm (7.2"); 10mm open



**SKU:** 240-071-03

**Categories:** [Open Surgery](#), [Neurosurgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Laminectomy Punches</a>
<b>Aesculap</b>	FF720
<b>KLS Martin</b>	24-993-03
<b>AUTHORS</b>	<a href="#">KERRISON</a>
<b>Medicon</b>	57.08.43
<b>SPECIFICATIONS</b>	<a href="#">40° up cut</a> , <a href="#">without ejector</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Laminectomy Punches</a>
<b>Aesculap</b>	FF720
<b>KLS Martin</b>	24-993-03

## AUTHORS

[KERRISON](#)

## Medicon

57.08.43

## SPECIFICATIONS

[40° up cut](#), [without ejector](#)