

KERRISON PUNCH; 40° UP CUT

3mm; without ejector; Kerrison ring hand

18.0cm (7.2"); 10mm open




SKU: 240-371-03

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
AUTHORS	KERRISON
SPECIFICATIONS	40° up cut , without ejector



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
AUTHORS	KERRISON
SPECIFICATIONS	40° up cut , without ejector