

KERRISON PUNCH; 40° UP CUT

5mm; without ejector; Kerrison ring hand

18.0cm (7.2"); 10mm open




SKU: 240-371-05

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
AUTHORS	KERRISON
SPECIFICATIONS	40° up cut , without ejector



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
AUTHORS	KERRISON
SPECIFICATIONS	40° up cut , without ejector