

KERRISON PUNCH; 40° UP CUT

4mm; with ejector; thin footplate

18.0cm (7.2"); 10mm open; CC-coated




SKU: 240-071-04ATCC

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
AUTHORS	CLEANWAVE , KERRISON
SPECIFICATIONS	40° up cut , with ejector



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
AUTHORS	CLEANWAVE , KERRISON
SPECIFICATIONS	40° up cut , with ejector